



Product Code	35936
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Yes
CH-REP	Yes

EAN/UPC	8023279359367
---------	---------------

Description

EURONDA EUROSEAL® INFINITY SEALING MACHINE

Euroseal® Infinity is a practical, sturdy, reliable and easy to use electronic thermosealing machine for continuous cycle processing at constant temperature.

Compact design

Compact and functional design to minimise the space required for its use so it can be comfortably placed on a table top or wall-mounted (need optional kit) in a vertical position.

Versatile

The roll holder is adjustable in two positions: high for best performance when the machine is wall-mounted, or low to speed up roll selection and changing operations.

Ergonomic

The new built-in cutter and sealing lever positioned at the front, not only improve ergonomics and safety, but also make packaging and sealing operations easier.

Safe

It has a metal bar that holds and locks the roll and ensures that it remains inside the sealing area. The smooth surfaces make it easy to clean.

Reliability

The quality of the materials and sturdy design make it indestructible, guaranteeing long-lasting efficiency and reliability.

Compliant with UNI EN 868-5 and EN 11607-2.

Multilingual manual: GB, IT.

Available QR code manual in: FR, ES, PT, DE, PL, CZ, SI, RO, RS, RU.

Technical Specifications

Sealing size: 12x310 mm
Sealing temperature: 170-200°C
Size - weight: 430x430xh 215 mm - 7.1 kg
Power supply: 200-240 V; 50/60 Hz, 100W