



Product Code	35939
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	
UK-REP	Yes
CH-REP	Yes
EAN/UPC	8023279359398

Description

EURONDA EUROSEAL® SEALING MACHINE

Euroseal® is a compact, ergonomic and easy to use thermosealing machine.
The see-through polycarbonate cover allows to perform sealing operations with greater precision.

Versatile
Its special patented roll holder is adjustable in height and depth and can be wall-mounted (need optional kit) in order to meet the specific needs of each work space.

Practical
The integrated support surface aids the sealing process and provides more spacious and practical working conditions.

Safe
Anti roll-back blocking system and safety system to prevent the roll from burning.
Audible and visual warning signals during use.

Sturdy
Lightweight and cutting-edge materials, built to withstand heavy usage.

Fully compliant with UNI EN 868-5 and EN ISO 11607-2.
Multilingual manual: GB, IT.
Available QR code manual in: FR, ES, PT, DE, CZ, SE.

Technical Specifications

Sealing size: 12x310 mm
Sealing temperature: 170-200°C
Size - weight: 456x390xh 315/420* mm - 7.1 kg
Power supply: 200-240 V; 50/60 Hz, 100W

*When roll holder is extended (double roll holder)