



Product Code	35952
Unit of sale	box of 1000 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279359527

Description

TYPE 5 MOVING FRONT INTEGRATOR

- box of 1,000 pcs
- non-toxic / lead free

Moving front integrators 134 °C - 3.5 min are giving you the best interpretation possible in case of any failure in your steam sterilization cycle.

If temperature is not reached or holding time was too short the dark colour will not reach the accept area.

Conform to ISO 11140 - type 5.

User manual: GB, IT, FR, ES.

- shelf life: 5 years after production
- packed in white aluminium re.closable bags
- very easy interpretation of result
- complete indicator area sealed of by foil