



Product Code	35953
Unit of sale	box of 500 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279359534

## Description

### TYPE 5 STEAM INTEGRATOR

- box of 500 pcs
- non-toxic / lead free

Integrators can be used in addition to the batch control Helix which allow to release load.  
It allow to re-affirm sterilization condition inside trays or containers.  
SV: 134 °C - 3.5 min / 121 °C - 15 min.  
Conform to ISO 11140 - type 5.

Package and manual in: GB, IT, FR, ES, PT, GR, DE.

- shelf life: 5 years after production
- packed in an easy dispensing carton
- clear colour transition from blue into pink of lacquer covered indicator