



Product Code	35957
Unit of sale	box of 250 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC 8023279359572

Description

SEAL CHECKER

- box of 250 pcs
- dimension: 175 x 75
- suitable for all rotary heat sealers

The ISO-11607 standard is specifying a routine, daily check of your heat sealers. The seal checker is a way to perform the test required and to document which sealer, which day, which person did perform the test and to easily show results of your heat sealer capacity of properly heat seal your pouch. When place inside a pouch and seal, the sealed area will be darker and can be easily read whether the seal is consistent, regular and completely without any gaps.