



Product Code	35962
Unit of sale	box of 28 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	0707387583965

Description

3M™ Comply™ Lead free indicator tape 1322 - 18 mm x 55 m - box of 28 rolls

It consists of adhesive, backing and chemical indicator stripes.

Aggressive and pressure-sensitive adhesive designed to adhere to a variety of wraps, including untreated woven and disposable non-woven, paper and plastic wraps in order to secure the pack during steam sterilization.

Backing is a beige crepe paper that provides the stretch needed for pack expansion.

Chemical indicator lines shows a visual color change from off-white/tan to dark-brown/black when exposed to sterilization.

- not made with Natural Rubber Latex
- trusted, high quality adhesive
- lead Free Ink means no hazardous waste issues
- stretchable backing that minimizes tape "pop-off" during sterilization
- can be written on or labeled with pre-printed labels