

3M ATTEST INSTRUMENT PROTECTOR 5x12.5 cm - 13911



| | |
|---------------|-------------------|
| Product Code | 35964 |
| Unit of sale | box of 1000 pcs |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|------|----------------|
| GTIN | 50707387792003 |
|------|----------------|

Description

3M™ Attest™ instrument protector 5x12.5 cm 13911 - box of 1000

Designed to cover the pointed and cutting edges of surgical instruments in steam or ethylene sterilization. Reduces the possibility of pouch puncture during sterilization. Indicators on the protectors verify steam or EO gas processing Lead level.

Scissors not included.