

2 LINE STEAM LABEL for labeller 35980



Product Code	35984
Unit of sale	box of 12000 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8034049914648
PARAF	964968046

Description

2 LINE STEAM LABEL - BOX OF 12,000 (24 ROLLS OF 500 LABELS)

First print line can display: 3 alphabetic digits + 11 numbers Second print line can display: 14 numbers Steam indicator label, suitable for steam sterilization cycle.

Double adhesive support. After the sterilization cycle the indicator changes from original PINK colour to BROWN.

Process indicator: Steam Standard ISO 11140-1 Class 1, non toxic ink, latex free.