

STEAM ADHESIVE LABEL 32x19mm for labeller 35985 - DM Class I



Product Code	35986
Unit of sale	box of 12000 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Not required
CH-REP	Not required

RDM (NSIS)	2322848
CND	Z1399
EAN/UPC	8034049916154
PARAF	990689186

Description

<div>2 LINE STEAM LABELS - box of 12,000 (24 rolls of 500 labels)</div><div></div><div>Labels with double adhesive support and steam indicator. When exposed to steam sterilization, the process indicator turns from pink to brown.</div><div>Label size: 32x19 mm.</div><div>Steam indicator: steam standard ISO 11140-1 type 1.</div><div>No latex and no heavy metal.</div><div></div>