



| | |
|---------------|-------------------|
| Product Code | 35987 |
| Unit of sale | box of 5 pcs |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

UK-REP Not required

CH-REP Not required

EAN/UPC 8034049916116

Description

<div>INK ROLL FOR LABELLER - box of 5</div><div></div><div>Tamp/Print H ink is ideal for medical labeller.</div><div></div><div>Packed in blister of 5 pieces.</div><div>Nickel free.</div><div>It complies with EN ISO 14001:2015 e EN ISO 9001:2015</div><div></div>