


COMPATIBLE PAEDIATRIC PADS for defibrillator
Medtronic Physio Control see also 55004



Product Code	33603
Unit of sale	kit of 2 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2469016
CND	C020401
EAN/UPC	8033695728265
GMDN	47055

Description

Pediatric Compatible pads - disposable
(CPR guidelines state 8 years as limit age for pediatric treatment, corresponding to about 25 kg according to standard growth curves)

These pads have an energy reduction system: defibrillator output of 200J is reduced to 50J.
International compliance to regulation as above adult compatible pads.
Latex free.
Made in Italy.

Multilingual envelope: GB, FR, IT, ES, PT, DE, NL, SE, GR, RU.

Compatible with:

- Physiocontrol
- Lifepak 500 AED pink socket
- Lifepak 1000 AED
- Lifepak CR plus AED
- Lifepak Express AED