

SELF-SEAL POUCHES 70x230 mm



Product Code	36031
Unit of sale	box of 4000 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes
RDM (NSIS)	2180053
CND	S01010101
GTIN	18033844701115
GMDN	13735

Description

SELF-SEAL POUCHES FOR STERILIZATION BY STEAM OR ETO

External size: 70 x 230 mm
Minimum order: 4,000 pcs (Packed in box of 200 pouches)

Self-seal pouches provide very strong and equally divided sealings (multiline).
The pouches have an adhesive band with the innovative system “finger-lift” which allows an easier use.
Duo-check system: double sterilization indicator, internal and external.
It ensures high control in the sterilization process.

Indicators provide clear colour changes when exposed to sterilization process:

- steam: light blue to green
- ethylene oxide: pink to orange
- formaldehyde: green to yellow.