

PICCOLO PLUS DIAGNOSTIC AUDIOMETER - air+bone+mask+Maestro software



Product Code	33615
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2351779
CND	Z121401
EAN/UPC	8054187380372

Description

Piccolo Plus diagnostic audiometer - air + bone + mask
for Pc use

Piccolo is a basic diagnostic audiometer offering air and bone audiometry functions. Piccolo can be controlled from a computer running Windows. Lightweight, compact design, powerful, userfriendly software and Noah connectivity combine to make Piccolo the ideal choice for professionals on the move.

- computer control through the Maestro software.

Compatible for Windows 7, 8, 10

- small and light weight

- power supply through the USB connection, when controlled by computer.

Software language: GB, FR, IT, ES, DE, PL, PT, RU.

CD manual and quick guide GB, FR, IT, ES, DE, PL, PT, NL.

Piccolo needs annual calibration by authorized centre. Please check on www.gimaitaly.com list of countries where you can find a centre.

Made in Italy.

Technical Specifications

Type of device: Pc control
Class EN 60645-1, ANSI S3.6: 3-diagnostic
Output: AC, BC
Masking: NBN
Warble tone: Yes
Pulsed tone: 0.5 Hz, 1 Hz, 2 Hz
Frequency range: 125 - 8,000 Hz
Maximum output intensities: AC: -10 to 120 dB HL / BC: -10 to 80 dB HL
Attenuator step: 5 dB
Automatic audiometry: HW auto threshold
Talk over: External microphone
Connection to PC: USB [driverless]
Controlling software: Maestro module
Controlling app: - Maestro (App Store)
Dimensions: 16x16x3cm (16.3x6.3x1.2 in)
Weight: 300 g (10.6 oz)
Operating voltage: 100-240 V, 50/60 Hz

Standard accessories

- DD445 supra-aural headphones
- B71 bone vibrator
- Patient response switch
- Clip-on microphone for operator to-patient communication
- Inventis Daisy Patient Record Manager, with Daisy Maestro module
- Noah Maestro module
- USB connection cable
- Medical grande power supply
- User Manual