



Product Code	33619
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2351781
CND	Z121401
EAN/UPC	8054187380778

Description

TRIANGLE BASIC SCREENING AUDIOMETER - air

Triangle is a battery-operated, portable screening audiometer, capable of performing automatic and manual pure tone audiometry exams in air conduction. Designed to deliver fast, accurate, and reliable pure tone audiometry for healthcare professionals on the move. Its unique combination of advanced technology and userfriendly design, featuring a capacitive colour touchscreen, allow for intuitive operation, making hearing screening simpler and more efficient than ever.

- Pure Tone Audiometry in AC
- compact and lightweight
- rechargeable lithium battery delivers over 12 hours of continuous use
- internal memory for 100 patients
- Made in Italy.

Technical Specifications

Type of device: screening audiometer
 Class IEC 60645-1/ANSI S3.6: type 4
 Output: AC
 Warble tone: Yes
 Pulsed tone: 1 Hz, 2 Hz
 Frequency range: 125-8,000 Hz
 Maximum output intensities: AC: -10 to 100 dB HL
 Attenuator step: 5 dB
 Automatic audiometry: available algorithms: Hughson Westlake, Quick Search, Fixed Intensity
 Talk over: built-in microphone
 Compatible software: Maestro
 Memory: up to 100 patients
 Display: TFT LCD touchscreen display, 240x320 resolution
 Display size: 2.8" (73.8x57.6 mm)
 Power supply: built-in lithium rechargeable battery
 Operating voltage: 100-240 V, 50/60 Hz
 Size -weight: 21x16x4 cm - 425 g

Standard accessories

DD45 supra-aural headphones
 USB connection cable
 Patient response switch
 Soft carrying case
 Audiogram cards
 Medical grade power supply
 User manual: GB, FR, IT, ES, DE.
 On-line manuals
 PT, NL, DK, SE, NO, LT, LV, EE, PL, HU, RO, BG, GR, TR available on gimaitaly.com