

# ETHICON VICRYL PLUS ABSORBABLE SUTURES - gauge 2/0 needle 24 mm - braided



Product Code	22363
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	2487113
CND	H010101040201
EAN/UPC	8023279223637
GTIN	30705031078264
GMDN	45401

## Description

VICRYL PLUS BRAIDED - absorbable sterile surgical sutures

Complete absorption in 56-70 days

- Ethicon code: VCP453H
- Gauge: 2/0
- Metric: 3.0
- Colour: violet
- Circle: 3/8
- Needle length: 24 mm
- Needle: FS-1
- Suture length: 70 cm
- Minimum order: box of 36 pieces

vicryl plus is indicated for use in general tissue approximation and/or ligation except for ophthalmic, neuro, cardiovascular.