

# ANTIMICROBIAL ADHESIVE FILM 105 Mic - roll 5m x 400 mm



Product Code	36505
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279365054
PARAF	986898423

## Description

ANTIMICROBIAL adhesive film 105 mic - roll 5 m

- SILVER FILMS - roll 5 m x 400 mm
- Type: Adhesive
- Thickness: 105 micron
- Colour: Clear (Anti-Glare)
- Durability: > 6 months
- Applications: Touchscreens, IT devices, displays, ...
- Medical Device Directive 93/42/ECC
- Antiviral test: ISO21702 Antibacterial test JIS Z-2801

Silver film developed to reduce indirect transmissions of virus and bacteria, available in three versions: clear adhesive and non-adhesive, silver coloured adhesive.

It reduces the viral load of 99,24% within 30 minutes, 99,76% within 1 hour and 99,97% within 24 hours, after coming into contact with virus and bacteria.

Each type can be applied on a wide range of objects and surfaces in different environments: public and medical facilities, home, shops, devices.

Unique 105 mic or 200 mic thickness ensures more efficacy compared to similar products, the 105 mic or 200 mic thickness and silver presence provide better protection.

Silver-ion prevents bacteria and virus proliferation: silver nanoparticles penetrate inside the virus destroying bacteria through a response called "IONIC reaction"

High transparency and anti-glare

Anti-fingerprints and anti-shatter: special coating eases removal of fingerprints or stains and protects from damages.

Entirely made of PET which maintains clear and clean status.

3H level of Hardness assures durability (>6 months) and protection from external shocks and scratches.