



Product Code	36542
Unit of sale	box of 30 pcs
Minimum order	1
Type	Medical device
Class	III

RDM (NSIS)	2374137
CND	A02010701
EAN/UPC	8023279365429

Description

BD POSIFLUSH XS PREFILLED SYRINGE 10 ml - 306572 - box of 30 - sterile

BD PosiFlush™ XS Prefilled Saline Syringe is a sterile, ready to use, polypropylene syringe with Luer-lock™ tip, that contains sterile and non-pyrogenic isotonic 0.9% sodium chloride solution. It is used in oncology, interventional radiology and critical care for flushing of in-situ peripheral intravenous catheters (PIVCs), peripherally inserted central catheters (PICCs), central venous catheters (CVCs) and implanted venous access parts.

Reduction of occlusions and catheter damage

Designed to reduce the likelihood of syringe-induced catheter reflux, that can cause occlusions. 10 ml diameter helps to lower the risk of catheter damage due to injection pressure

IV catheter-related infections reduction

Luer-lock™ tip cap helps to prevent touch contamination and infection and retaining ring reduces the risk of solution contamination

Clear labelling helps to reduce the risk of medication errors

Easily identifiable syringe contents and bold print for clarity

No need for needles, glass vials and ampoules

It can be used on a sterile field by healthcare professionals to reduce the risk of touch contamination that may occur during manual preparation of saline flush syringes.

It is not formulated with natural rubber latex. Class III

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, NO, FI, DK, SE, GR.