

B-BRAUN ECOTAINER STERILE IRRIGATION SOLUTION - 500 ml



Product Code	36555
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	395138
CND	V9099
PARAF	977815671
GTIN	04030539077200

Description

B-BRAUN ECOTAINER® - 500 ml

Suggested for quick irrigation with high volume of solution.

Multilanguage manual and label: GB, FR, IT, ES, DE, PT, NL, DK, FI, NO, SE.

Composition (1,000 ml solution contain):

- sodium chloride 9 g
- water
- Injection electrolytes: mmol/l
- Na+ 154
- Cl- 154