

B-BRAUN OMNIFLUSH SYRINGE - 10 ml



Product Code	36558
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	483561
CND	V9099
EAN/UPC	4046964043215
PARAF	977815683

Description

B-BRAUN OMNIFLUSH® - 10 ml

PREFILLED FLUSH syringe 10 ml with NaCl 0.9% - box of 100 - EM-3513576

The Omniflush® syringe is designed specifically to minimize unintended blood reflux into the catheter lumen thus reduces the risk of catheter blockage. With Omniflush® the fluid amount used for flushing can be reduced to the required minimum. For flushing of all intravenous accesses.

Multilanguage manual and label: GB, FR, IT, ES, DE, PT, PL, CZ, BG, DK, EE, FI, GR, HU, LT, LV, NL, NO RO, SK, SI, SE.

Composition (1,000 ml solution contain):

- sodium chloride 9 g
- water
- Injection electrolytes: mmol/l
- Na+ 154
- Cl- 154