

POVI-IODINE 100 ANTISEPTIC - 125 ml - biocide - marketable in IT only



Product Code	36560
Unit of sale	box of 24 pcs
Minimum order	1
Type	Biocide
Class	

UK-REP Yes

EAN/UPC 8034105280045

Description

POVI-IODINE 100 antiseptic - 125 ml - box of 24

Povi Iodine 100 is a ready to use liquid solution based on 10% concentration of Iodopovidone. Ideal for hand washing and disinfection in healthcare facilities.

Professional use.

Biocide effective action

Bactericide (EN 13727, EN 1500, EN 1499), Yeasticidal (EN 13624), Virucide (Adenovirus, Murine noravirus EN 14476)

Biocide PT01 - IT/2021/00714/aut.

Registered and marketable in Italy only.

Label: IT.

Shelf life: 24 months.

iodo-POVIDONE

10% IODINE

Made in Italy.

Technical Specifications

Chemical composition

100 g of solution contain: 10 g Polyvinylpyrrolidone iodine (iodo -povidone 10% iodine), glycerin, coformulants and purified water q.b. 100 g