

# POVI-IODINE 100 ANTISEPTIC - 500 ml - biocide - marketable in IT only



|               |               |
|---------------|---------------|
| Product Code  | 36565         |
| Unit of sale  | 1 pc          |
| Minimum order | 1             |
| Type          | Biocide       |
| Class         |               |
| UK-REP        | Yes           |
| EAN/UPC       | 8034105280069 |

## Description

POVI-IODINE 100 antiseptic - 500 ml

Povi Iodine 100 is a ready to use liquid solution based on 10% concentration of Iodopovidone. Ideal for hand washing and disinfection in healthcare facilities.

Professional use.

Biocide effective action

Bactericide (EN 13727, EN 1500, EN 1499), Yeasticidal (EN 13624), Virucide (Adenovirus, Murine noravirus EN 14476)

Biocide PT01 - IT/2021/00714/aut.

Registered and marketable in Italy only.

Label: IT.

Shelf life: 24 months.

ODO-POVIDONE  
10% IODINE

Made in Italy.

## Technical Specifications

Chemical composition

100 g of solution contain: 10 g Polyvinylpyrrolidone iodine (iodo -povidone 10% iodine), glycerin, coformulants and purified water q.b. 100 g