

GERMOCID MULTIUSO SPRAY DISINFECTANT DETERGENT - 750 ml



Product Code	36579
--------------	-------

Unit of sale	box of 12 pcs
--------------	------------------

Minimum order	1
---------------	---

Type	PMC - disinfectant
------	--------------------

Class	
-------	--

UK-REP	No
--------	----

CH-REP	No
--------	----

EAN/UPC	8009110030671
---------	---------------

GTIN	08009110030671
------	----------------

Description

GERMOCID MULTIUSO - 750 ml - box of 12 bottles

Germocid Multiuso is a ready to use disinfectant detergent, specially designed for cleaning delicate surfaces. It can be used in the bathroom and kitchen, on tiles and all washable surfaces. It is active against Gram+ and Gram- bacteria and fungi when left for 15 minutes. It degreases, deodorizes and cleanses ensuring a complete disinfection of objects and environments. It dries quickly and leaves no streaks. It does not foam and stain.

Multilingual label: GB, FR, IT.
Italian MOH Registration No 19.259.

Made in Italy.

Technical Specifications

Chemical composition

100,00 g of product contain: Benzalkonium Chloride 0,50 g; O-phenilphenol 0,02 g; coformulants; depurated water q.s. to 100,00 g.