



Product Code	36600
Unit of sale	box of 600 pcs
Minimum order	1
Type	PMC - disinfectant
Class	

PARAF	986899375
GTIN	18009110030357

Description

TEKNA disinfectant wipes for surfaces - box of 600 (12 flowpack of 50)

Ready to use wipes designed to cleanse and disinfect simultaneously objects and surfaces, eliminating germs and bacteria. The presence in the formulation of Chlorhexidine, Benzalkonium and Alcohol guarantees the deep disinfection of all washable surfaces, with an excellent detergent and degreasing action. Do not contain parfum, colorants and whitening agents.

In accordance with EN1650, EN1276, EN13697. CE regulation N 648/2004.

Label: GB, IT.

Wipe size: 14.3x23.3 cm

Made in Italy.

Technical Specifications

Chemical composition

100.0 g of imbibed solution of alcoholic disinfectant contains: Chlorhexidine Digluconate 0.1 g; Benzalkonium chloride 0.1 g; Isopropyl alcohol; Ethyl alcohol; coformulants; depurated water q.s. to 100.0 g. Italian MOH Registration No. 19354