



x 12

Product Code	36605
Unit of sale	box of 12 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
PARAF	987469071

Description

OXYGENATED WATER (Hydrogen peroxide) - 1 l - box of 12

Hydrogen peroxide (3% diluted solution) 10 volumes.

Clear solution, colourless.

Multilanguage label: GB, FR, IT, ES.

Made in Italy.

Technical Specifications

CHEMICAL COMPOSITION:

Hydrogen peroxide: =3% (10 vol.)

Water: =96.999%

Stabilizer: 0.001%