

GERMOCID BASIC 750 ml without vaporizer



Product Code	36610
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	4926
CND	D99
GTIN	08009110029781

Description

GERMOCID BASIC - 750 ml - box of 12 (without vaporizer)

Alcohol-based disinfectant cleaner for non-invasive medical devices.

Disinfectant action: the association between Chlorhexidine and Quaternary ammonium Salt and the high alcoholic content guarantees a large action spectrum towards:

- Gram+ and Gram- bacteria, fungi (yeasts and molds) and viruses (HIV, HBV, HCV included) in 5 minutes.
- mycobacteria (TBC included) in 60 minutes.

Strong detergent power and fast evaporation. It is not de-activated by organic substances. Bactericide, virucide and fungicide activity.

Label/instruction: GB, FR, IT

Made in Italy.

Technical Specifications

Chemical composition

100 g of product contain: Chlorhexidine digluconato 0,1 g; Benzalkonium Chloride 0,1 g; Ethyl Alcohol 45,0 g; Isopropyl Alcohol 15,0 g as needed to 100 g

Product form: Colourless liquid. No added perfume.