



Product Code	36616
Unit of sale	1 pc
Minimum order	1
Type	PMC - disinfectant
Class	

EAN/UPC	8009110028319
---------	---------------

## Description

### ENVIRONMENT DISINFECTANT - BOTTLE OF 1 LITRE

Concentrated liquid detergent-disinfectant for environments. Dilute in water to a ratio of 0.5% for the disinfection and cleaning of surfaces, counters, tiles, floors, and sanitary fixtures, including when using floor-washing machines. Suitable for large areas such as hospitals, clinics, rest homes, and gyms.

- Disinfectant action: wide protection against viruses, fungi, and Gram+ and Gram- bacteria.
- Detergent action: contains non-ionic surfactants to enhance surface cleaning.

Label/instruction: GB, FR, IT.

Made in Italy.

## Technical Specifications

Chemical composition:

100 g OF PRODUCT CONTAIN:

Quaternary ammonium compound: 10 g

Isothiazolinone magnesium chloride: 0.2 g

Coformulants and deperated water: 100 g

APPEARANCE: Pale yellow liquid

SMELL: Citrus

BIODEGRADABILITY: Detergent according to Regulation CE N° 648/2004

ITALIAN HEALTH MINISTRY AUTHORIZATION N°: 9658