



Product Code	36618
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I
RDM (NSIS)	2588281
CND	D0801
EAN/UPC	8009110026193

## Description

### ENZYME-BASED DETERGENT - 3 I

Medical device enzymatic based detergent for invasive surgical instruments.

This clear liquid detergent is very useful to improve the cleaning of orthodontic and surgical instruments, glassware, plastic and metal objects. Also suitable for ultrasound cleaners.

Dilute 0.6-10% in 1 litre warm water, immerse the instruments for 10 minutes and rinse.

Label/instruction: GB, FR, IT.

Made in Italy.

## Technical Specifications

Chemical composition (100 g contain): proteolytic enzymes 0.5 g, non-ionic surfactant agents 10-15 g, coformulants, water as needed to 100 g