



Product Code	36621
Unit of sale	box of 12 pcs
Minimum order	1
Type	
Class	

EAN/UPC	8009110030275
PARAF	986465603

Description

Office spray disinfectant - 150 ml - box of 12

Ready to use disinfectant spray suitable for quick and effective disinfection of the surfaces and equipment in the workplaces: keyboards, phones, desks, calculators and more. Acts in 15 minutes.

It does not stain, does not wet, does not damage tissues and materials and releases a pleasant scent.

Active on the Coronavirus family.

Tested by an accredited EU laboratory.

Bactericidal activity: EN 1040; EN1276; NF T 72-281 EN 13697; EN 13727; Fungicidal activity: EN 1275; EN 1650; EN 13697

Virucidal activity (HIV, HBV, HCV, Coronavirus): EN14476.

Label: GB, IT

Made in Italy.

Technical Specifications

Chemical composition

100,00 g of the product contains: Benzalkonium chloride 0,10 g; Orthophenylphenol 0,02 g; Isopropyl alcohol; parfum; propellant q.s. to 100,00 g.

Italian MOH Registration No. 7604