



Product Code	36623
Unit of sale	1 pc
Minimum order	1
Type	PMC - disinfectant
Class	

EAN/UPC	8009110050020
PARAF	987476482

## Description

### TEKNA DISINFECTANT SPRAY - 400 ml

Disinfectant spray approved by Italian Health Ministry.

In addition to quickly eliminating fungi, bacteria and allergens, it performs an effective deodorant action.

Useful to disinfect and deodorize soft fabrics and surfaces such as sofas, armchairs, chairs, cushions, mattresses, rugs, curtains, carpets or high passage areas such as dressing rooms, waiting rooms, classrooms, canteens, gyms, public toilets.

#### Disinfectant action:

Complete action spectrum against Gram-positive and Gram-negative bacteria and fungi.

Tekna disinfectant spray does not wet, does not stain, does not leave residues and does not damage fabrics and materials.

## Technical Specifications

### CHEMICAL COMPOSITION

100,00 g contain: Benzalkonium chloride 0,10 g; Biphenyl-2-ol 0,02 g; Propan-2-ol; Perfume; Propellants as needed to 100,00 g.

Product form: perfumed spray. Italian MOH Registration No. 13.431