



Product Code	36629
Unit of sale	box of 12 pcs
Minimum order	1
Type	PMC - disinfectant
Class	

EAN/UPC	8023279366297
PARAF	987469828

Description

Medical SOAP - 0.5 l - box of 12

Antiseptic liquid for hands and skin disinfection. It is indicated for surgeon and medical staff's hands and also patient skin disinfection. It reduces risk of cross-infections. It contains emollient substance, which make skin soft.

Label/instructions GB, IT.

Made in Italy.

Technical Specifications

Chemical composition

100.0 g contain: Alkyl benzyl dimethyl ammonium chloride 0.50 g; O-phenylphenol 0.50 g; perfume; coformulants, deperated water q.s. to 100.0 g.

Italian Health Ministry Registration 19570