



| | |
|---------------|--------------------|
| Product Code | 36631 |
| Unit of sale | box of 16 pcs |
| Minimum order | 1 |
| Type | PMC - disinfectant |
| Class | |

| | |
|--------|----------------|
| UK-REP | No |
| CH-REP | No |
| GTIN | 08009110029972 |

Description

FLOWPACK 20 WIPES WITH CHLORHEXIDINE - display of 16 flowpacks

Disinfectant TNT wipes with chlorhexidine are skin antiseptic, indicated for disinfection of hands and intact skin in general. Fast skin disinfection scented with bergamot oil. Ready-to-use with broad-spectrum bacterial action.

Label/instructions GB, IT

Italian MOH Registration No. 11.571

Made in Italy.

Technical Specifications

CHEMICAL COMPOSITIONS

100 g of solution contain: Chlorhexidine digluconate 1 g, soluble bergamot essential oil, deperated water q.s. to 100 g, limonene, linalool