

# MULTI USE DISINFECTANT WIPES WITH CHLOREXIDINE



Product Code 36633

Unit of sale 1 pc

Minimum order 1

Type PMC - disinfectant

Class

UK-REP No

CH-REP No

EAN/UPC 8009110028906

## Description

DISINFECTANT WIPES WITH CHLORHEXIDINE - tube of 100 wipes

Disinfectant TNT wipes with chlorhexidine are skin antiseptic, indicated for disinfection of hands and intact skin in general. Fast skin disinfection scented with bergamot oil. Ready-to-use with broad-spectrum bacterial action.

Label/instructions GB, IT

Italian MOH Registration No. 11.571

Made in Italy.

## Technical Specifications

### CHEMICAL COMPOSITIONS

100 g of solution contain: Chlorhexidine digluconate 1 g, soluble bergamot essential oil, depurated water q.s. to 100 g, limonene, linalool