



Product Code	36665
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A
RDM (NSIS)	1873521
CND	D99
EAN/UPC	8009110026230
PARAF	987475151

Description

GERMOCID inodore - tank 5 l

Concentrated solution with a rapid and effective action as detergent and disinfectant.

It is active against Gram+ and Gram- bacteria, fungi, mycobacteria and viruses (Adenovirus, HIV, HBV, HCV).

Rapid and practical cleansing and disinfection of all surfaces of non-invasive medical devices such as beds, litters, masks and various equipment (desks, armrests, handles, etc.) usually in use at medical or dental surgeries, hospitals, beauty centers.

It disinfects, degreases and cleanses, does not stain and does not leave halos.

Also suitable for alcohol-sensitive surfaces and materials.

Class IIa Medical Device in compliance with Directive 93/42/EEC.

Label/instruction: GB, IT.

Made in Italy.

Technical Specifications

Chemical composition

100,0 g of product contain:

Benzalkonium Chloride 10,0 g;

Orthophenylphenol 1,0 g;

Coformulants, depurated water q.s. to 100,0 g