

STAND FOR HAND SANITIZER DISPENSER AND GLOVES



Accessori non inclusi

Accessories not included

Accessoires non inclus

Product Code 36676

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279366761

Description

STAND FOR HAND SANITIZER DISPENSER AND GLOVE BOX

Sturdy and resistant stand made of epoxy painted steel, black colour.
It can be positioned in any type of room.
Base disassembled with feet.

Made in Italy.

Sign holder

Holder for glove box (gloves not included)

Holder (base 10x10 cm) for hand sanitizer dispenser (dispenser not included)

Base Ø 60 cm