



Product Code	36677
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8009110026148

Description

DOSASOAP NO-TOUCH HAND WASH SYSTEM - white

Wall-mounted no-touch dispenser with photocell for antiseptics, convenient for detergents, liquid soaps and gels. A sensor detects hand movements and dispenses the product. It ensures a rapid and safe hand hygiene, without any contact with the dispenser.

Suitable for 1 litre bottles (36630, 36632, 36644, 36595, 36671, 36672).

Battery operated, it needs 4xAA, not included.