



Product Code	36681
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	1872848
CND	D99
EAN/UPC	8009110019508

## Description

GERMOCID TEC - 750 ml

Foaming disinfectant cleaner for non-invasive medical devices used in tattoo studios, clinics, hospitals, and beauty shops. Suitable for sensitive surfaces where alcohol-based disinfectants cannot be used. It is a safe product that dries quickly, does not stain, and leaves no streaks.

Ready-to-use spray that cleans and, thanks to the synergic action of the Quaternary Ammonium Salt with Orthophenylphenol, disinfects against bacteria, yeasts, fungi, moulds, TBC, and viruses (including Adenovirus, HIV, HBV, and HCV).

Label/instruction: GB, FR, IT.

Made in Italy.

## Technical Specifications

### CHEMICAL COMPOSITION.

100 g of product contain: Benzalkonium Chloride; O-phenylphenol; coformulants; deperated water as needed to 100 g;

Product form: Colourless liquid with citron perfume