


# DECONTAMINATING MAT 90x115 cm - 30 layers - blue



Product Code	36695
Unit of sale	box of 5 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2449893
CND	T030501
EAN/UPC	8023279366952
UMDNS	11740
PARAF	977815986
GTIN	08023279366952
GMDN	44011

## Description

Mats to be used in all environments that require total absence of contaminating agents and the preservation of high hygienic levels.

GIMA mat, thanks to antibacterial glue captures and retains dust and bacteria and substantially reduces bacteria growth and spreading.

A full length adhesive sheet on the bottom of the mat provides a firm to grip to the floor. Easy removability without leaving any adhesive trace.

Each mat includes 30 numbered polythene sheets and non-adhesive edge strip to allow easy removal of each used layer.

Packed in a stylish dispenser folder with handle for easy carriage and vertical storage. Each mat can be easily picked leaving the others perfectly stored.

- Total mat thickness: 1.5 mm
- Dimensions: 90 x 115 cm
- Color: blue
- Box of 5 mats with 30 sheet each