

SOFFISOF AIR DRY INCONTINENCE PAD - moderate incontinence - large



Product Code	36708
Unit of sale	box of 90 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2414304
CND	T04010103
EAN/UPC	8007300011172
PARAF	972068112
GTIN	28007300009609

Description

SOFFISOF AIR DRY

- Size: Large
- Circumference (cm): 100-150
- Moderate incontinence-day
- Minimum order: box of 90 (6 boxes of 15)

Anatomically shaped absorbent diapers with repositionable fastening tapes.
Made of fluff of pure cellulose, superabsorbent polymers and waterproof exterior layer.
They ensure maximum safety in case of moderate or high incontinence level.
Hypoallergenic, dermo protective with odour control.
The new design guarantees maximum comfort.
Gradual fade wetness indicator easily suggests when it's time for a change.
Air Dry pads have a special breathable outer layer that helps to reduce the onset of redness and irritations.

Multilingual package: GB, FR, IT, ES, DE, PT, Arabic.
Made in Italy.