

SOFFISOF AIR DRY INCONTINENCE PAD - heavy incontinence - XXL



Product Code	36712
Unit of sale	box of 60 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2414428
CND	T04010103
EAN/UPC	8007300011233
PARAF	972968187
GTIN	28007300009685
GMDN	11239

Description

SOFFISOF AIR DRY

- Size: XX-Large
- Circumference (cm): 155-175
- Heavy incontinence-night
- Minimum order: box of 60 (4 boxes of 15)

Anatomically shaped absorbent diapers with repositionable fastening tapes.
Made of fluff of pure cellulose, superabsorbent polymers and waterproof exterior layer.
They ensure maximum safety in case of moderate or high incontinence level.
Hypoallergenic, dermo protective with odour control.
The new design guarantees maximum comfort.
Gradual fade wetness indicator easily suggests when it's time for a change.
Air Dry pads have a special breathable outer layer that helps to reduce the onset of redness and irritations.

Multilingual package: GB, FR, IT, ES, DE, PT, Arabic.
Made in Italy.