



Product Code	36719
Unit of sale	box of 32 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2358861
CND	T04010103
EAN/UPC	8007300028675
PARAF	977610474
GTIN	28007300028679
GMDN	11239

Description

SOFFISOF AIR DRY PANTS/PULL-UP

- Size: Large
- Circumference (cm): 110-150
- Heavy incontinence-night
- Minimum order: box of 32 (4 boxes of 8)

Disposable absorbent panty ideal for walking, self-sufficient and active patients, with slight to severe leakages.

Air dry, hypoallergenic, dermo protective with odour control.

Outer lining in soft breathable fabric allows skin to breathe.

Ergonomic elastic band for maximum comfort and greater fit.

Soft elastic and thigh line barriers ensure maximum leak protection.

Gradual fade wetness indicator easily suggests panty change.

Multilingual package: GB, FR, IT, ES, DE, PT, Arabic.

Made in Italy.