

# PREMIUM PATIENT ANTI-HANDLING SUIT - XL - reusable



Product Code	36733
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279367331

## Description

Patient anti-handling suit

- Size: XL
- Unisex
- Green

The suit prevents the patient from handling himself during night or day.  
The zip allows the operator to intervene easily on the patient, for hygiene.  
Back closure with double slider zip and final reinforcement.  
Long sleeves with soft elastic at the cuffs and bottom of the leg; double reinforced crotch Made of 100% polyester honeycomb stretch weft.  
Resistant to washing with chlorine at 38°C and in water with detergent at 90°C.  
Weight: gr/m2 200 ±5%.  
Quick drying and thermoregulating.

>100 WASHES: Resistant to washing with chlorine at 38°C and in water with detergent at 90°C

Made in Italy.