

PATIENT ANTI-HANDLING SUIT - L - reusable



Product Code	36738
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279367386

Description

PATIENT ANTI-HANDLING SUIT

Unisex - green
Size: L
Height x chest: h 154x128 cm

The suit prevents the patient from handling Back himself during night or day. Back closure with double slider zipper and final reinforcement. The zip allows the operator to intervene easily on the patient, for hygiene. Long sleeve and leg with elastic cuffs and ankles; double reinforced crotch. Designed for elderly people with incontinence problems associated with neurodegenerative diseases, clothing intolerance or mental disabilities. Made of jersey: 97% cotton and 3% lycra. Washable up to 60°C. Weight: 160 gr/m2.