



Product Code	36801
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
--------	-----

RDM (NSIS)	2226019
CND	D0799
EAN/UPC	4028163096209

Description

INCIDIN™ LIQUID - tank 5 l

Ready to use surface disinfectant based on alcohol.

Ideal for:

- disinfection of medical devices surfaces
- intensive care and high risk areas. It offers high user safety thanks to its broad spectrum of effectiveness which also includes TB and Noro.
- acts against Bacteria, Yeasts and Tuberculosis in 5 minutes
- reduced contact time
- formulation free of dyes and perfume
- excellent material compatibility on alcohol resistant surfaces.

Medical Device, Class IIa.

Multilingual label: GB, FR, IT, DE, NL, CZ, SK, HR, RS.

Technical Specifications

Chemical composition.

100 g contain: 35g 2-propanolo, 25g 1-propanolo