



Product Code	36805
Unit of sale	box of 6 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
RDM (NSIS)	2200918
CND	D0599
EAN/UPC	4028163115832

Description

INCIDIN™ OXYFOAM S - spray 750 ml - box of 6

Ready to use sporicidal detergent and disinfectant with Hi-Speed H2O2™

Ideal for:

- cleaning and disinfection of medical devices
- terminal disinfection of semi-critical medical devices (probes, TEE probes, ...)

The advanced formulation provides the benefits of a patented solution of hydrogen peroxide with an effective broad spectrum.

Effective against Bacteria, Yeasts, Fungi and Mycobacteria, avoiding active residues thanks to the rapid evaporation rate.

OxyFoam S provides:

- active ingredients decomposing in water and oxygen
- sporicidal and virucidal efficacy
- excellent material compatibility

Medical Device, Class IIb

Multilingual label: GB, FR, IT, DE, NL, CZ, SK, HR, RS, EE, LT, LV.

Technical Specifications

Chemical composition

100 g contain: 1.5 g hydrogen peroxide