



Product Code	36810
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	1441611
CND	D99
EAN/UPC	8009110025509

Description

GERMOCID INSTRUMENTS - DILUITED 2% - bottle 1 l - box of 12

Ready to use solution with a detergent and disinfectant action for invasive medical devices and surgical instruments. Specific product for the cleaning, highlevel disinfection and decontamination before the sterilization of surgical instruments and invasive medical devices. It has a high efficiency as detergent and disinfectant also in presence of organic material (blood, saliva, organic exudates) and it can be also used on dental instruments. The dirt removal does not require mechanical action. It is active against Gram+ and Gram-bacteria, fungi, mycobacteria and viruses (HIV-HBV-HCV).

For professional use.

Medical Device class IIb in compliance with Directive 93/42/EEC.

Instructions: GB, FR, IT.

Made in Italy.

Technical Specifications

CHEMICAL COMPOSITION

100 g of product contain: Benzalkonium chloride 2,0 g; O-phenylphenol 0,2 g; coformulants; depurated water q.s. to 100,0 g.