

# CONTAINER WITH FILTER small h 135 mm - grey



Product Code	37001
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2811854
CND	S9099
EAN/UPC	8023279370010

## Description

Standard container with one filter and non perforated bottom.

- Small model: 285 x 280 mm
- Height: 135 mm
- Non perforated bottom
- 1 filter
- Colour: grey

Standard containers (supplied with 1 paper filter and filter holder).  
All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.