

CONTAINER WITH FILTER medium h 100 mm - red



Product Code	37060
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC	8023279370607
---------	---------------

Description

Standard container with one filter and non perforated bottom.

- Medium model: 465 x 280 mm
- Height: 100 mm
- Non perforated bottom
- 1 filter
- Colour: red

Standard containers (supplied with 1 paper filter and filter holder).
All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.