

# CONTAINER WITH FILTER medium h 135 mm - blue - perforated



Product Code	37081
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC 8023279370812

## Description

Standard container with two filters and perforated bottom.

Perforated bottom allows much quicker sterilization process (half the time compared to not perforated). Lock nail on bottom side made of chrome.

Standard containers (supplied with 1 paper filter and filter holder)

- Medium model: 465 x 280 mm
- Height: 135 mm
- Perforated bottom
- 2 filters
- Colour: blue

All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.