

# CONTAINER WITH FILTER medium h 150 mm - blue - perforated



Product Code	37082
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2989805
CND	S9099
EAN/UPC	8023279370829

## Description

Standard container with two filters and perforated bottom.

Perforated bottom allows much quicker sterilization process (half the time compared to not perforated). Lock nail on bottom side made of chrome.

Standard containers (supplied with 1 paper filter and filter holder)

- Medium model: 465 x 280 mm
- Height: 150 mm
- Perforated bottom
- 2 filters
- Colour: blue

All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.