

ETHICON VICRYL RAPID ABSORBABLE SUTURES - gauge 4/0 needle 19 mm - braided



Product Code	22371
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2436964
CND	H0101010299
EAN/UPC	8023279223712
GTIN	30705031106547
GMDN	17471

Description

VICRYL RAPID BRAIDED - absorbable sterile surgical sutures

Complete absorption in 42 days

- Ethicon code: V2920G
- Gauge: 4/0
- Metric: 1.5
- Colour: undyed
- Circle: 3/8
- Needle length: 19 mm
- Needle: FS-2
- Suture length: 45 cm
- Minimum order: box of 12 pieces

Vicryl Rapid Braided indicated only for: use in superficial soft tissue approximation of the skin and mucosa, where only short term (7-10 days) wound support is required.